

Specifications**SoundSleep CPAP/AutoCPAP**

Mode	LC-CPAP LC-APAP	CPAP CPAP, Auto
Pressure Range		4 to 20 cmH ₂ O (hPa)
Pressure Accuracy		± (0.5 hPa+4%)
Ramp Time		0 to 60 Min (5-min increments)
Starting Ramp Pressure		4 cmH ₂ O to CPAP Pressure (CPAP); 4 cmH ₂ O to Min Pressure (Auto CPAP)
Expiratory Pressure Relief		0 to 3
Leak Compensation		Automatic
Altitude Compensation		Automatic
Humidification		0 to 5
Data Storage Capacity		SD card & On-board: > 1 year
Dimensions		14x15.6x9.4cm (w/o humidifier) 26.2x15.6x9.4cm (w/ humidifier)
Weight		950g (w/o humidifier) 1500g (w/ humidifier)
Electrical Requirements		100-240 VAC, 50/60 Hz, 2.0- 1.0 A
Noise Level		<30 dB @ 10 cmH ₂ O
Data Analysis Software		PC Based
Filter		Reusable Pollen
Tubing		22mm, 6ft. (1.8m)
Device Controls		Control Dial/Push Button
Language		Chinese/English/Spanish
Operating Temperature		5° to 35° C (41° to 95° F)
Storage Temperature		-20° to 60° C (-4° to 140° F)
Relative Humidity (operating & storage)		10 to 93% (non-condensing)
Atmospheric Pressure		101 to 77 kPa (0-2286 m/0-7500 ft)
Risk Classification		Class IIa, Type BF, IP22