

cuffs

Overlapping chambers system

CarePump devices work with 4- / 6- / 8-chamber cuffs, designed for lymphatic drainage of legs, arms, hips, and the area covering the abdomen, bottom and thighs.

Features

- + overlapping system (the chambers in the cuff overlap, which prevents fluid retreat during the treatment)
- + zipper along the entire cuff
- + heavy duty tubes
- + high-quality material, easy to clean
- + cuffs have an universal size - after filling with the air they perfectly adapt to the body



leg cuffs



arm, shoulder and chest cuff



short pants cuff



arm cuff



waist cuff



full pants cuff

comparison of models

	EXPERT8	MOVE8 ^{PRO}	COMPACT4	LITE8	LITE6	LITE4	MOVE8	MOVE6	MOVE4
Line/series	Advanced	Advanced	Advanced	Economic	Economic	Economic	Sport	Sport	Sport
Target group	B2B	B2B	B2B/B2C	B2C HomeCare	B2C HomeCare	B2C HomeCare	B2C Sport	B2C Sport	B2C Sport
Number of chambers	8	8	4	8	6	4	8	6	4
Max. pressure [mmHg]	220	220	220	220	220	220	220	220	220
Therapeutic pre-programs	13	15	6	8	6	4	6	6	6
Certificate of the Polish Lymphological Society	✓	✓	✓	✓	✓	✓			
Graphical presentation of pre-programs progress	✓	✓	✓	✓	✓	✓	✓	✓	✓
User's own programs	✓	✓	✓						
List of favourite programs	✓	✓	✓						
Actual pressure measurement in the chambers	✓	✓	✓	✓	✓	✓	✓	✓	✓
Precise regulation of pressure, time and speed in steps of 1, 5, 10	✓	✓	✓	✓	✓	✓	✓	✓	✓
Setting individual pressure in each chamber	✓	✓	✓				✓	✓	✓
'Smart Gradient' function	✓	✓	✓				✓	✓	✓
'Interval' function	✓	✓	✓						
'Hold' function	✓	✓	✓						
Bag included		✓					✓	✓	✓
Battery powered		✓					✓	✓	✓
Features for visually impaired people	✓		✓						
Change of language	✓	✓	✓	✓	✓	✓	✓	✓	✓
Voice commands	✓		✓						
Display size [inch]	7	7	5	5	5	5	7	7	7
Weight [kg]	3,6	2,9	3,1	3,2	3,1	2,9	2,9	2,8	2,5
Dimensions [cm]	30x26x16	28,5x12,5x12,5	30x26x16	30x26x16	30x26x16	30x26x16	28,5x12,5x12,5	28,5x12,5x12,5	28,5x12,5x12,5

CARE PUMP

BardoMed Sp. z o.o.
2F/6 Pod Fortem Str.
31-302 Krakow, Poland

www.carepump.com
/CarePumpPL
@carepump_official

CARE PUMP



compression therapy equipment
lymphatic drainage

advanced line

Advanced interface and the highest functionality



- + 8- (Expert8, Move8PRO) / 4- (Compact4) chamber cuffs
- + Max. pressure: 220 mmHg
- + 15 (Move8PRO) / 13 (Expert8) / 6 (Compact4) therapeutic pre-programs of the medical organization
- + Target group: B2B/B2C

- + User's own programs
- + List of 'Favourite' programs
- + Modification of 'Ready and quick save' program
- + Smart Gradient' function (mmHg, %, final pressure)
- + Adjustment of 'Interval' function – a pause between cycles
- + Adjustment of 'Hold' function – keeping pressure in the chambers
- + Variable interval for setting the parameter value
- + Setting individual pressure in each chamber
- + Reverse pumping option

- + Ergonomic navigation keyboard
- + Features for visually impaired / blind people: voice assistant + special keyboard
- + Change of language and voice commands with adjustment / off option
- + Sound of the device being ready to work/end of therapy with the adjustment / off option
- + Extensive settings and service menu
- + Password protection on startup of the device
- + Time code protection of the device's working time
- + Patient safety button option

sport line

Recovery packed to go



- + 8- (Move8) / 6- (Move6) / 4- (Move4) chamber device
- + Max. pressure: 220 mmHg
- + 6 therapeutic pre-programs
- + Target group: B2C Sport

- + Smart Gradient' function (mmHg, %, final pressure)
- + Setting individual pressure in each chamber
- + Reverse pumping option
- + Dedicated backpack included
- + Battery-powered

MOVE8^{PRO}

the most advanced battery-powered device

The perfect combination of the highest functionality of Expert8 device with portable design from CarePump Sport series. 3-hours battery life allows to make a treatment anywhere, with the most customized parameters.

15 message modes



economic line

Simplicity and functionality



- + 8- (Lite8) / 6- (Lite6) / 4- (Lite4) chamber device
- + Max. pressure: 220 mmHg
- + 8 (Lite8) / 6 (Lite6) / 4 (Lite4) therapeutic pre-programs of the medical organization
- + Target group: B2C HomeCare
- + Knob for parameters adjustment
- + Patient safety button option

All CarePump devices features

- + Certificate of a medical organization
- + Actual pressure measurement in the chambers
- + Precise regulation of pressure, time and speed in steps of 1
- + Change of parameters during treatment
- + Adjustment of the speed of the treatment / massage
- + Real time view of the details and therapy progress
- + Graphical presentation of pre-programs progress
- + Disable / bypass selected chambers
- + Modern design
- + Glass panel with an integrated color touch screen