

CARE PUMP




LITE4 *Economic line*

Number of chambers	4
Max. pressure [mmHg]	220
Number of proprietary therapeutic pre-programs of the medical organization	4
Precise regulation of pressure and time	in steps of 1
Number of independent channels	1
Adjustment of the speed of the treatment	1-5
Cuffs	overlapping system
Display	5" TN lcd color touch screen
Weight [kg]	2,9
Dimensions [mm]	300x260x160
Power supply	100÷240 V, 50÷60 Hz
Power consumption	40-70 W
Warranty	1 year

Accessories




Included:




power supply

Optional:




4-chamber arm cuff

CPARM4T




4-chamber waist cuff

CPWAIST4T



4-chamber left leg cuff

CPLEG4TL




4-chamber right leg cuff

CPLEG4TR




4-chamber short pants

CPSHORTPANTS4T




Extension zipper

CPLEGEXT4




Double splitter

CPDOUBLECONNECTOR4




Patient's safety button

SAFETY BUTTON



Double splitter (long - 150 cm)
for 4-chamber devices

CPDOUBLECONNECTOR4L



2 x Extension zippers
for 4-chamber short pants

CPSHORTPANTSEXT4

Therapeutic pre-programs

of the medical organization



lymphoedema I

Time	50 min.
Hold	3 s
Interval	3 s

lymphoedema II

Time	50 min.
Hold	3 s
Interval	3 s

prophylaxis, relaxation

Time	60 min.
Hold	2 s
Interval	2 s

all-chamber massage

Time	30 min.
Hold	2 s
Interval	3 s