

# CARE PUMP




## **LITE6** *Economic line*

|                                                                            |                              |
|----------------------------------------------------------------------------|------------------------------|
| Number of chambers                                                         | 6                            |
| Max. pressure [mmHg]                                                       | 220                          |
| Number of proprietary therapeutic pre-programs of the medical organization | 6                            |
| Precise regulation of pressure and time                                    | in steps of 1                |
| Number of independent channels                                             | 1                            |
| Adjustment of the speed of the treatment                                   | 1-5                          |
| Cuffs                                                                      | system overlapping           |
| Display                                                                    | 5" TN lcd color touch screen |
| Weight [kg]                                                                | 3,1                          |
| Dimensions [mm]                                                            | 300x260x160                  |
| Power supply                                                               | 100÷240 V, 50÷60 Hz          |
| Power consumption                                                          | 40-70 W                      |
| Warranty                                                                   | 1 year                       |

# Accessories




## Included:




power supply

## Optional:




6-chamber arm cuff

CPARM6T




6-chamber waist cuff

CPWAIST6T



6-chamber left leg cuff

CPLEG6TL




6-chamber right leg cuff

CPLEG6TR



6-chamber short pants

CPSHORTPANTS6T




Extension zipper

CPLEEXT6




Double splitter

CPDOUBLECONNECTOR6




Patient's safety button

SAFETY BUTTON




6-chamber double cuff for arms, shoulders and chest

CPDOUBLEARM6T




6-chamber left arm, shoulder  
and chest cuff

CPARMSB6TL




6-chamber right arm, shoulder  
and chest cuff

CPARMSB6TR




6-chamber full pants

CPLEG6TL




2 x Extension zippers  
for 6-chamber double cuff for  
arms, shoulders and chest

CPDOUBLEARMEXT6




2 x Extension zippers  
for 6-chamber short pants

CPSHORTPANTSEXT6



2 x Extension zippers  
for 6-chamber full pants

CPLONGPANTSEXT6



Double splitter (long - 150 cm)  
for 6-chamber devices

CPDOUBLECONNECTOR6L

# Therapeutic pre-programs

## of the medical organization



lymphoedema I

- 
- 
- 
- 
- 
- 

|          |         |
|----------|---------|
| Time     | 50 min. |
| Hold     | 3 s     |
| Interval | 3 s     |

lymphoedema II

- 
- 
- 
- 
- 
- 

|          |         |
|----------|---------|
| Time     | 50 min. |
| Hold     | 3 s     |
| Interval | 3 s     |

lipoedema I

- 
- 
- 
- 
- 
- 

|          |         |
|----------|---------|
| Time     | 60 min. |
| Hold     | 3 s     |
| Interval | 3 s     |

regeneration

- 
- 
- 
- 
- 
- 

|          |         |
|----------|---------|
| Time     | 45 min. |
| Hold     | 2 s     |
| Interval | 2 s     |

prevention, relaxation

- 
- 
- 
- 
- 
- 

|          |         |
|----------|---------|
| Time     | 60 min. |
| Hold     | 2 s     |
| Interval | 2 s     |

all-chamber massage

- 
- 
- 
- 
- 
- 

|          |         |
|----------|---------|
| Time     | 30 min. |
| Hold     | 2 s     |
| Interval | 3 s     |