

CARE PUMP



LITE8


Economic line

Number of chambers	8
Max. pressure [mmHg]	220
Number of proprietary therapeutic pre-programs of the medical organization	8
Precise regulation of pressure and time	in steps of 1
Number of independent channels	1
Adjustment of the speed of the treatment	1-5
Cuffs	system overlapping
Display	5" TN lcd color touch screen
Weight [kg]	3,2
Dimensions [mm]	300x260x160
Power supply	100÷240 V, 50÷60 Hz
Power consumption	40-70 W
Warranty	1 year

Accessories




Included:




power supply

Optional:




8-chamber arm cuff

CPARM8T




8-chamber waist cuff

CPWAIST8T




8-chamber left leg cuff

CPLEG8TL




8-chamber right leg cuff

CPLEG8TR




8-chamber short pants

CPSHORTPANTS8T




Extension zipper

CPLEGEXT8




Double splitter

CPDOUBLECONNECTOR8




Patient's safety button

SAFETY BUTTON




8-chamber double cuff for arms, shoulders and chest

CPDOUBLEARM8T




8-chamber left arm, shoulder
and chest cuff

CPARMSB8TL




8-chamber right arm, shoulder
and chest cuff

CPARMSB8TR




8-chamber full pants

CPLONGPANTS8T




2 x Extension zippers
for 8-chamber double cuff for
arms, shoulders and chest

CPDOUBLEARMEXT8




2 x Extension zippers
for 8-chamber short pants

CPSHORTPANTSEXT8



2 x Extension zippers
for 8-chamber full pants

CPLONGPANTSEXT8



Double splitter (long - 150 cm)
for 8-chamber devices

CPDOUBLECONNECTOR8L

Therapeutic pre-programs

of the medical organization



mode **A**

-
-
-
-
-
-
-
-

Time	50 min.
Hold	3 s
Interval	3 s

mode **B**

-
-
-
-
-
-
-
-

Time	50 min.
Hold	3 s
Interval	3 s

mode **C**

-
-
-
-
-
-
-
-

Time	50 min.
Hold	3 s
Interval	0 s

mode **D**

-
-
-
-
-
-
-
-

Time	60 min.
Hold	3 s
Interval	3 s

mode **E**

-
-
-
-
-
-
-
-
-
-
-
-

x3

Time	60 min.
Hold	1 s
Interval	2 s

mode **F**

-
-
-
-
-
-
-
-

Time	45 min.
Hold	2 s
Interval	2 s

mode **G**

-
-
-
-
-
-
-
-

Time	60 min.
Hold	2 s
Interval	2 s

mode **H**

-
-
-
-
-
-
-
-

Time	30 min.
Hold	2 s
Interval	3 s