

What makes them stand out?



6 operating modes dedicated to athletes
(including innovative **mixed mode** enabling personalization of the procedure)



PROGRADIENT System – pressure drop in subsequent pumped chambers



SUPERFLOW technology – 5 speed levels



Real-time pressure measurement in the cuff chambers (range 20-220 mmHg)



Precise parameter adjustment



Possibility to disable any cuff chambers, e.g. to avoid injured areas



Reverse pumping system



Powerful and **efficient pump** with a lifespan of up to 10,000 h



Replaceable battery



A wide selection of cuffs for various body parts



Multilingual menu



Made in the EU, Polish origin

CARE PUMP

(+48) 722 000 017

export@carepump.com

www.carepump.com

 CarePump

 @carepump_official

 /CarePumpPL

CARE PUMP

Lymphatic drainage in sports

Effective post-workout recovery



Lymphatic drainage – support in regeneration

People practicing sports often struggle with **muscle fatigue, swelling** and **injuries** that make it difficult to return to regular training. Lymphatic drainage is one of the most effective methods of **post-workout regeneration, training support** and **injury prevention**.

How does lymphatic drainage work?

This is active regeneration combined with relaxation. The treatment involves placing a cuff on a given part of the body, the chambers of which, thanks to the connection to the pressotherapy device, are successively filled with air and, through sequential **compression**, produce a **massage effect**.

The treatment:

- + stimulates blood and lymph circulation
- + nourishes and oxygenates tissues
- + supports the removal of metabolic by-products, residual fluids and toxins
- + massages and relaxes muscles
- + accelerates the reconstruction of overloaded tissues

Benefits?

- + Reduction of muscle pain, tension and fatigue after exhausting training or competitions
- + Faster regeneration rate
- + Feeling of „light legs”
- + Increased body performance
- + Improvement of sports results
- + Less risk of injury
- + Reduction of oedema
- + Relaxation

Lymphatic drainage can be performed even daily. It is perfect for training camps and at the peak of the sports season, when the athlete needs to maintain **top form for a long time**.

CarePump **SPORT** series

CarePump **Sport** are lymphatic drainage devices designed for athletes, trainers and sports physiotherapists.

These are **mobile** and battery-powered devices distinguished by their modernity and extremely high effectiveness. They are compact and light, equipped with a convenient **backpack**, thanks to which athletes and their physiotherapists can easily take them to training camps or competitions. **Battery power** allows you to perform the treatment anywhere, also outdoors, immediately after training or sports competitions.

Available in 4-, 6- and 8-chamber versions.

The more chambers in the cuffs, the more accurate the massage and the greater the effectiveness.



MOVE4

4-chamber cuffs



MOVE6

6-chamber cuffs



MOVE8

8-chamber cuffs