

comparison of models

	EXPERT8	MOVE8 ^{PRO}	COMPACT4	LITE8	LITE6	LITE4	MOVE8	MOVE6	MOVE4
Line/series	Advanced	Advanced	Advanced	Economic	Economic	Economic	Sport	Sport	Sport
Target group	B2B	B2B	B2B/B2C	B2C HomeCare	B2C HomeCare	B2C HomeCare	B2C Sport	B2C Sport	B2C Sport
Number of chambers	8	8	4	8	6	4	8	6	4
Max. pressure [mmHg]	220	220	220	220	220	220	220	220	220
Therapeutic pre-programs	13	15	6	8	6	4	6	6	6
Certificate of the Polish Lymphological Society	✓	✓	✓	✓	✓	✓			
Graphical presentation of pre-programs progress	✓	✓	✓	✓	✓	✓	✓	✓	✓
User's own programs	✓	✓	✓						
List of favourite programs	✓	✓	✓						
Actual pressure measurement in the chambers	✓	✓	✓	✓	✓	✓	✓	✓	✓
Precise regulation of pressure, time and speed in steps of 1, 5, 10	✓	✓	✓	✓	✓	✓	✓	✓	✓
Setting individual pressure in each chamber	✓	✓	✓				✓	✓	✓
'Smart Gradient' function	✓	✓	✓				✓	✓	✓
'Interval' function	✓	✓	✓						
'Hold' function	✓	✓	✓						
Bag included		✓					✓	✓	✓
Battery powered		✓					✓	✓	✓
Features for visually impaired people	✓		✓						
Change of language	✓	✓	✓	✓	✓	✓	✓	✓	✓
Voice commands	✓		✓						
Display size [inch]	7	7	5	5	5	5	7	7	7
Weight [kg]	3,6	2,9	3,1	3,2	3,1	2,9	2,9	2,8	2,5
Dimensions [cm]	30x26x16	28,5x12,5x12,5	30x26x16	30x26x16	30x26x16	30x26x16	28,5x12,5x12,5	28,5x12,5x12,5	28,5x12,5x12,5